




ACUTE SURGICAL WOUND ANIMAL STUDY

HELIOGEN[®] vs. Type I Bovine Collagen Particulate^{*}



Objective: Evaluate the host response to scaffold products^{1,2}

Methods: Subcutaneous implant of scaffold products in athymic nude mice^{**} with acute wounds

Test Articles: Scaffold Products

HELIOGEN

- Type I and III bovine dermis particulate
- n=8 to 10 animals per time point[^]

Type I Bovine Collagen

- Bovine dermis particulate
- n=6 to 7 animals per time point[^]

Timepoints Evaluated:

- 1 week and 4 weeks
- Tissue sample section processed and evaluated microscopically


^{*} Not hydrolyzed collagen.

^{**} Athymic nude mice lack T cells. Therefore, they do not reject foreign tissue. This model eliminates confounding results from a xenogenic response.

[^] One site per animal.

Results: Week 1 Post-Implantation¹


HELIOGEN is **well-tolerated** as evident by **fewer inflammatory macrophages** infiltrating into the implant



Results: Week 4 Post-Implantation²

New collagen was deposited into the **HELIOGEN** implant site because the scaffold was **well-tolerated and supported cellular ingrowth**.

Type I Bovine Collagen Particulate was **largely degraded**; therefore, minimal collagen deposition was apparent within the implant.



Staining Guide



Macrophages (brown staining) are inflammatory cells responsible for clearing damaged tissue and foreign objects³



All nucleated cells (blue/purple staining)

¹Images are representative of observed response across all animals. Additional details are available upon request.

Conclusion^{1,2}

Week 1: HELIOGEN exhibited a lesser inflammatory response than the Type I Bovine Collagen Particulate.

Week 4: HELIOGEN persisted and exhibited scaffold properties like supporting new collagen production. The Type I Bovine Collagen Particulate was degraded due to early inflammation.

Overall, HELIOGEN collagen matrix provided an optimal scaffold in an acute surgical wound.[†]



Product	Product Compatibility With Host [‡]	Product Persistence Post-Implantation
HELIOGEN	++++	++++
Type I Bovine Collagen Particulate	++	++

+Ratings are based upon histopathological analysis. Higher ratings (+) are indicative of improved performance.



SCAN TO WATCH VIDEO PRESENTATION

LEARN MORE

To find out more about MIMEDX products:

Please Call: 866.477.4219 Email: customerservice@mimedx.com

[†]Results were demonstrated in animal only. Correlation to clinical results in humans is unknown.
[‡]Biocompatibility was determined by the macrophage response because nude mice lack and adaptive immune system.

1. MIMEDX Internal Report. MM-RD-00183, Immunofluorescence and Immunohistochemical Analysis of Xenograft ECM Particulates in a Mouse Subcutaneous Implant Model. 2. MIMEDX Internal Report. MM-RD-00178, Non-GLP Evaluation of Xenograft Scaffolds for Cellular Response in a Mouse Subcutaneous Implant Model. 3. Aitcheson SM, Frentiu FD, Hurn SE, Edwards K, Murray RZ. Skin Wound Healing: Normal Macrophage Function and Macrophage Dysfunction in Diabetic Wounds. Molecules. 2021;26(16):4917.

Note: Mouse and stain guide images created using BioRender.
Refer to Instructions for Use for a complete description of all applicable indications, warnings, precautions and contraindications.

