

# Evaluation of Lyophilized Human Amnion/Chorion Membrane (LHACM) in the management of nonhealing Diabetic Foot Ulcers (DFU)

An interim analysis of the CAMPAIGN Prospective Randomized Controlled Trial (RCT)

**BACKGROUND:** Chronic DFUs remain a significant challenge and are associated with high rates of infection, hospitalization, and amputation. The CAMPAIGN RCT was designed to assess whether LHACM, when added to a standardized Standard of Care (SOC) regimen, improves clinical outcomes.



N=88 randomized patients



12 Weeks



47% Closure Rates with LHACM plus SOC @ 12 Weeks



21% Closure Rates with SOC alone @ 12 Weeks

## STUDY DESIGN & METHODOLOGY

**Design:** Prospective, multicenter, RCT

**4-Week Run-In Period:** Subjects received standardized SOC during screening. Ulcers with  $\geq 20\%$  area reduction during run-in were excluded to ensure enrollment of non-healing DFUs

**Interim Analysis:** N=88, randomized patients

**Planned Enrollment:** N=150, randomized patients


**Comparison Arms:** LHACM + Standard of Care (SOC) compared to SOC alone. SOC Included:

- Sharp/mechanical debridement
- Moisture-appropriate dressing (calcium alginate or foam)
- Infection evaluation and management
- Mandatory offloading, including use of the Foot Defender® where appropriate

**Primary Endpoint:** Probability of complete closure at 12 weeks.

**Methodological Approach:** The CAMPAIGN trial used a Bayesian Framework which is explicitly supported by CMS.

## RESULTS



### At 12 weeks:

- 47% of wounds achieved complete closure with LHACM plus SOC
- 21% of wounds achieved complete closure with SOC

**Bayesian Framework:** LHACM exceeded the interim superiority threshold of  $>90\%$  with a posterior probability of 98.5%, supporting a high confidence of its clinical benefit.

**CONCLUSION:** The CAMPAIGN interim RCT demonstrated that LHACM provided a clinically meaningful advantage over SOC alone, nearly doubling closure probability by 12 weeks. These results reinforce the role of evidence-based amnion/chorion allografts in DFU management for high-quality clinical evidence.

To find out more about MIMEDX products:

Call: 866.477.4219 or email: [customerservice@mimedx.com](mailto:customerservice@mimedx.com)